

## Anti-H3N2 antigen from influenza A virus [IF1A11] Bulk Size Ab00443-3.0-BT

This chimeric mouse antibody was made using the variable domain sequences of the original human Fab format, for improved compatibility with existing reagents, assays and techniques.

Isotype and Format: Mouse IgG2b, Kappa

**Clone Number:** IF1A11

Alternative Name(s) of Target: Hemagglutinin; HA UniProt Accession Number of Target Protein: P03437

Published Application(s): Blocking, WB, ELISA Published Species Reactivity: Influenza A

Immunogen: RNA from peripheral blood lymphocytes of a healthy volunteer who had recovered from an

Influenza A virus infection.

**Specificity:** The antibody binds to Influenza A H3N2 with a Kd of 3.6 nM.

**Application Notes:** The antibody binds to H3N2, which binds to sialic acid-containing receptor on the host cell surface leading to attachment and internalization of the virion. The antibody neutralizes the H3N2-receptor interaction.

**Antibody First Published in:** Shibuya et al. Identification of a human monoclonal Fab with neutralizing activity against H3N2 influenza A strain from a newly constructed human Fab library. Microbiology and Immunology 2008; 52(3):162-170 PMID:18402598

Note on publication: Describes the generation of a Fab library and a H3N2-neutralizing antibody.

## **Product Form**

**Size:** 1 mg Purified antibody in bulk size. **Purification:** Protein A affinity purified

Supplied In: PBS only.

**Storage Recommendation:** Store at 4°C for up to 3 months. Note, this antibody is provided without added preservatives, it is therefore recommed this antibody be handled under sterile conditions. For longer

storage, aliquot and store at -20°C.

**Concentration:** 1 mg/ml.

Important note – This product is for research use only. It is not intended for use in therapeutic or diagnostic

| © 2024 Absolute Antibody      | https://absoluteantibody.com/product/anti-h3n2-antigen-from-influenza-a |
|-------------------------------|---|
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
|                               |   |
| procedures for humans or anim | nals.   |
|                               |   |
|                               |   |