

## **Product Datasheet**

Sales Enquires: <a href="mailto:sales@absoluteantibody.com">sales@absoluteantibody.com</a>
Support Queries: <a href="mailto:support@absoluteantibody.com">support@absoluteantibody.com</a>

Tel: +44 (0) 1865 920810 Fax: +44 (0) 1865 920811 absoluteantibody.com

## Recombinant Human Free Lambda Light Chain

Cat No: Pr00116

## **Product Summary**

Description: Recombinant human free lambda light chain manufactured using AbAb's Recombinant

Platform Protein: Human Free Lambda Light Chain

Fc domain: n/a

Structure / Form: Disulfide-linked homodimer

Species: Human

Construct Design Note(s): Lambda light chain anti-Beta-galactosidase [13R4] antibody (c.f.

Ab00135) Host: HEK293

UniProt Accession Number: n/a

Alternative Description: Published Application(s):

Tested Applications(s):

Activity:

## **Product Form**

Purification: Affinity Matrix

Supplied in: PBS with preservative (0.02% Proclin 300)

Endotoxin: <1.0 EU/mg

Shipping: The product is shipped on blue ice. Upon receipt, store it immediately at the temperature recommended.


**Storage Recommendation:** Store at 4°C for up to 1 month. For longer term storage aliquot in small volumes and store at -20°C. Avoid repeated freeze-thaw cycles.

SDS PAGE Purity: >95%, as determined by SDS-PAGE and visualized by Coomassie Brilliant Blue.

Calculated Molecular weight (dimer): 45540 Da

Extinction coefficient: 76610 (calculation performed as described by Pace et al. (1995), PMID: 8563639).

Important note - This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals



Coomassie-stained SDS-PAGE gel of human recombinant free lambda light chain. Starting from the tissue culture supernatant, a single purification step was performed to enrich for recombinant free light-chains, both as monomer and dimer (at around 28 and 50 kDa, respectively).

procedures for humans or animals
Cat No: Pr00116 Recombinant Human Free Lambda Light Chain
Pr00111-1.9 under non-reducing and reducing (DTT) conditions resolved by SDS-PAGE and stained using Coomassie-Blue.
Important note - This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals

Important note - This product is for research use only. It is not intended for use in therapeutic or diagnostic